

Fructosamina

LOT

1906319680

LOT

Nivel 1

319000

LOT

Nivel 2

319010

VALORES ASIGNADOS - TECNICA MANUAL / ANALIZADORES						
TECNICA/ ANALIZADOR	NIVEL 1			NIVEL 2		
	VALOR MEDIO	RANGO		VALOR MEDIO	RANGO	
Konelab Series / Wiener lab, CT Series	237 umol/L	190	284	374 umol/L	299	449
	2.1 mmol/L (DMF)	1.6	2.5	3.3 mmol/L (DMF)	2.6	3.9
Targa - BT 3000	226 umol/L	181	271	392 umol/L	314	470
	2.0 mmol/L (DMF)	1.6	2.4	3.4 mmol/L (DMF)	2.7	4.1
Metrolab 2300 plus / Wiener lab. CM Series	220 umol/L	176	264	374 umol/L	299	449
	1.9 mmol/L (DMF)	1.5	2.3	3.3 mmol/L (DMF)	2.6	3.9
BT 3000 plus / Wiener lab. CB Series	226 umol/L	181	271	392 umol/L	314	470
	2.0 mmol/L (DMF)	1.6	2.4	3.4 mmol/L (DMF)	2.7	4.1
CMD 800 Series	237 umol/L	190	284	381 umol/L	305	457
	2.1 mmol/L (DMF)	1.6	2.5	3.3 mmol/L (DMF)	2.6	4.0
Técnica manual	239 umol/L	191	287	454 umol/L	363	545
	2.1 mmol/L (DMF)	1.7	2.5	4.0 mmol/L (DMF)	3.2	4.7

Glosario/Glossário/Glossary/Glosariusz

Español

Português

English

Polski

nivel
valores asignados
rango
analizadores

nível
valores estabelecidos
faixa
analisadores

level
assigned values
range
analyzers

poziom
zakres wartości
zakres
analizator



Wiener lab.