

Fructosamina

Control - 2 niveles

LOT 1908330600

LOT Nivel 1 329620

LOT Nivel 2 329610

VALORES ASIGNADOS - TECNICA MANUAL / ANALIZADORES						
TECNICA/ ANALIZADOR	NIVEL 1			NIVEL 2		
	VALOR MEDIO	RANGO		VALOR MEDIO	RANGO	
Konelab Series /	229 umol/L	183	274	483 umol/L	386	580
Wiener lab, CT Series	2.0 mmol/L (DMF)	1.6	2.4	4.2 mmol/L (DMF)	3.4	5.0
Targa - BT 3000	221 umol/L	177	266	504 umol/L	403	605
	1.9 mmol/L (DMF)	1.5	2.3	4.4 mmol/L (DMF)	3.5	5.3
Metrolab 2300 plus /	225 umol/L	180	270	502 umol/L	401	602
Wiener lab. CM Series	2.0 mmol/L (DMF)	1.6	2.4	4.4 mmol/L (DMF)	3.5	5.2
BT 3000 plus /	221 umol/L	177	266	504 umol/L	403	605
Wiener lab. CB Series	1.9 mmol/L (DMF)	1.5	2.3	4.4 mmol/L (DMF)	3.5	5.3
CMD 800 Series	214 umol/L	171	257	487 umol/L	389	584
	1.9 mmol/L (DMF)	1.5	2.2	4.2 mmol/L (DMF)	3.4	5.1
Técnica manual	220 umol/L	176	264	601 umol/L	480	721
	1.9 mmol/L (DMF)	1.5	2.3	5.2 mmol/L (DMF)	4.2	6.3

Glosario/Glossário/Glossary/Glosariusz

<u>Español</u>	<u>Português</u>	<u>English</u>	<u>Polski</u>
nivel	nível	level	poziom
valores asignados	valores estabelecidos	assigned values	zakres wartości
rango	faixa	range	zakres
analizadores	analizadores	analyzers	analizator

