

IgE Control

Turbitest AA

LOT 1907328450

LOT Nivel 1 328420

LOT Nivel 2 328430

VALORES ASIGNADOS

BT 3000 PLUS / WIENER LAB. CB SERIES						
KIT	NIVEL 1			NIVEL 2		
	VALOR MEDIO (UI/ml)	RANGO ACEPTABLE		VALOR MEDIO (UI/ml)	RANGO ACEPTABLE	
IgE Control	62.8	44.0	81.6	416	312	520

KONELAB SERIES / WIENER LAB. CT SERIES						
KIT	NIVEL 1			NIVEL 2		
	VALOR MEDIO (UI/ml)	RANGO ACEPTABLE		VALOR MEDIO (UI/ml)	RANGO ACEPTABLE	
IgE Control	63.5	47.6	79.4	397	298	496

METROLAB 2300 PLUS / WIENER LAB. CM SERIES						
KIT	NIVEL 1			NIVEL 2		
	VALOR MEDIO (UI/ml)	RANGO ACEPTABLE		VALOR MEDIO (UI/ml)	RANGO ACEPTABLE	
IgE Control	74.7	56.0	93.4	395	296	494

WIENER LAB. CMD SERIES						
KIT	NIVEL 1			NIVEL 2		
	VALOR MEDIO (UI/ml)	RANGO ACEPTABLE		VALOR MEDIO (UI/ml)	RANGO ACEPTABLE	
IgE Control	65.9	49.4	82.4	397	298	496

Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
kit	kit	kit	zestaw
nivel 1	nível 1	level 1	poziom 1
nivel 2	nível 2	level 2	poziom 2
valores asignad	valores estabelecidos	assigned values	zakres wartości
valor medio	valor médio	mean	średnia wartość
rango acceptable	faixa aceitável	acceptable range	pożądany zakres

