

# IgE Control

## Turbitest AA

LOT 2101411480

LOT Nivel 1 411050

LOT Nivel 2 411060

### VALORES ASIGNADOS

BT 3000 PLUS / WIENER LAB. CB SERIES						
KIT	NIVEL 1			NIVEL 2		
	VALOR MEDIO (UI/ml)	RANGO ACEPTABLE		VALOR MEDIO (UI/ml)	RANGO ACEPTABLE	
IgE Control	69.4	48.6	90.2	411	308	514

KONELAB SERIES / WIENER LAB. CT SERIES						
KIT	NIVEL 1			NIVEL 2		
	VALOR MEDIO (UI/ml)	RANGO ACEPTABLE		VALOR MEDIO (UI/ml)	RANGO ACEPTABLE	
IgE Control	64.3	48.2	80.4	410	308	513

METROLAB 2300 PLUS / WIENER LAB. CM SERIES						
KIT	NIVEL 1			NIVEL 2		
	VALOR MEDIO (UI/ml)	RANGO ACEPTABLE		VALOR MEDIO (UI/ml)	RANGO ACEPTABLE	
IgE Control	78.0	58.5	97.5	398	299	498

WIENER LAB. CMD SERIES						
KIT	NIVEL 1			NIVEL 2		
	VALOR MEDIO (UI/ml)	RANGO ACEPTABLE		VALOR MEDIO (UI/ml)	RANGO ACEPTABLE	
IgE Control	67.9	50.9	84.9	403	302	504

### Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
kit	kit	kit	zestaw
nivel 1	nível 1	level 1	poziom 1
nivel 2	nível 2	level 2	poziom 2
valores asignad	valores estabelecidos	assigned values	zakres wartości
valor medio	valor médio	mean	średnia wartość
rango acceptable	faixa aceitável	acceptable range	pożądany zakres

