

# Proti U/LCR

## Control - 2 niveles

LOT 2002360500

LOT Nivel 1 358070

LOT Nivel 2 358080

VALORES ASIGNADOS PARA TECNICA MANUAL Y ANALIZADORES						
TECNICA/ ANALIZADOR	NIVEL 1			NIVEL 2		
	VALOR MEDIO (mg/dL)	RANGO (mg/dL)		VALOR MEDIO (mg/dL)	RANGO (mg/dL)	
Konelab Series / Wiener lab, CT Series	30.2	22.6	37.7	100	80	120
Metrolab 2300 plus / Wiener lab. CM Series	31.8	23.9	39.8	96.5	77.2	115.8
Targa BT 3000	32.5	24.4	40.6	97.4	73.0	121.7
BT 3000 plus / Wiener lab. CB Series	32.5	24.4	40.6	97.4	73.0	121.7
Wiener lab. CMD Series	30.7	23.0	38.4	99.2	79.4	119.1
Técnica Manual	34.2	25.7	42.8	98.6	73.9	123.2

### Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
nivel	nível	level	poziom
fecha	data	date	data
valores asignados	valores estabelecidos	assigned values	zakres wartości
valor medio	valor médio	mean	średnia wartość
rango	faixa	range	zakres
técnica manual	técnica manual	manual technique	metody manualne
analizadores	analisadores	analyzers	analizator

